



## ALOKA SSD-5500 ULTRASCHALLGERÄT

**SKU: R122012**

**5.800,00 €**

Phased Array Sonde UST-5286-2.5Mhz

Aloka UST-5280-5 TEE Sonde

Printer P91 &

Videorekorder SVO-9500MDP

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall / Doppler / Sonden](#)  
[/Zubehör](#)

## **PRODUKTBESCHREIBUNG**

Phased Array Sonde UST-5286-2.5Mhz

Aloka UST-5280-5 TEE Sonde

Printer P91 &

Videorekorder SVO-9500MDP

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg



[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN

**Weight** 150 kg



## GALLERIE

