



**ALOKA UST-52108 // UST 52108 CARDIAC  
TRANSDUCER|SSD-5000|SSD-5500  
SV|ALPHA 5|ALPHA 6|ALPHA 7|ALPHA 10**

**SKU: CY082020-VK**

[Aloka Ust-52108 Phased Array Cardiac  
Transducer| SSD-5000 | SSD-5500 Sv | Alpha 5 |  
Alpha 6 |Alpha 7 | Alpha 10](#)

Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg

www.mtc-hamburg.com  
info@mtc-hamburg.com  
Mobile 0049 176 205 800 86  
Tel / Homezone 0049 40 18 20 79 03  
Tel / Büro 0049 40 3 61 61 - 256  
Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall / Doppler / Sonden  
/Zubehör](#)

## PRODUKTBESCHREIBUNG

Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg

www.mtc-hamburg.com  
info@mtc-hamburg.com  
Mobile 0049 176 205 800 86  
Tel / Homezone 0049 40 18 20 79 03  
Tel / Büro 0049 40 3 61 61 - 256  
Fax. 0049 40 3 61 61 - 257



## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

