



BAXTER IMMUCENT II ZENTRIFUGE MIT TIMER

SKU: B042014

399,00 €

Typ. B50565-3G

220v // 280 Watt // 50Hz // WK 48Nm // 2900U/Mn // Fm.

12 x 10 g

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com


info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Laborequipment](#)

PRODUKTBESCHREIBUNG

Typ. B50565-3G

220v // 280 Watt // 50Hz // WK 48Nm // 2900U/Mn // Fm. 12 x 10 g

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com



Mobile. 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

