



[CF-Q160 I CF-Q160L KOLOSKOP](#) [VIDEOKOLOSKOP](#)

SKU: AD042024

~~4.999,00 €~~ **4.444,00 €**

zzgl. [Versandkosten](#)

Zustand: Einsatzfähig – Sehr guter Zustand

Lagerort: Hamburg



Categories: [Flexible Endoskope & Zubehör](#)

PRODUKTBESCHREIBUNG

- Channel Width: 3.7
- Working Length: L: 1680
- Field of View: 140
- Direction of View: Forward viewing Depth of Field: 3-100 Outer Diameter: 12.8 Outer
- Diameter Insertion Tube: L: 12.8
- Max Angulation Up: 180
- Max Angulation Down: 180
- Max Angulation Right: 160
- Max Angulation Left: 160 Key Benefits

ZUSÄTZLICHE INFORMATIONEN

Weight

5 kg



GALLERIE

