



ERBE APC 300 ARGON BEAMER ARGON PLASMA KOAGULATION

SKU: AV102016 // EY112017-R1 /
A112018/DZ032019

777,00 €

zzgl. [Versandkosten](#)

Typ. 10132-010

Anwendung : in Kombination mit ERBOTOM ACC
450, 451, ICC 200 Argon Coag, 300, 350.

Anwendung : Zur Argon-Plasma-Koagulation

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [HF-Chirurgie & Zubehör](#)

PRODUKTBESCHREIBUNG

Typ. 10132-010

Anwendung : in Kombination mit ERBOTOM ACC 450, 451, ICC 200 Argon Coag, 300, 350.

Anwendung : Zur Argon-Plasma-Koagulation

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile

0049 176 205 800 86

Tel / Homezone

0049 40 18 20 79 03



Tel / Büro

0049 40 3 61 61 - 256

Fax.

0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN

GALLERIE

