



[FUJINON BRO-Y3S BRONCHOSCOPE](#)

SKU: S082015

800,00 €

zzgl. [Versandkosten](#)

Insertion Tube: 5.7 mm, Length: 57,5 cm, Total Length: 77 cm,

Channel Size: 2 mm

Field of View: 100°

Depth of Field: 3 - 50 mm

Angulation Up/Down: 160°/100°

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Endoskopie und Zubehör](#)

PRODUKTBESCHREIBUNG

Insertion Tube: 5.7 mm, Length: 57,5 cm, Total Length: 77 cm,

Channel Size: 2 mm



Field of View: 100°
Depth of Field: 3 - 50 mm
Angulation Up/Down: 160°/100°

Vollfunktionsfähig
Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com

Mobile. 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

