



## FUJINON EC-200 MR // EC 200 MR KOLOSKOP

**SKU:** BI102016-01

Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg

www.mtc-hamburg.com  
info@mtc-hamburg.com  
Mobile 0049 176 205 800 86  
Tel / Homezone 0049 40 18 20 79 03  
Tel / Büro 0049 40 3 61 61 - 256  
Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Flexible Endoskope & Zubehör](#)

## PRODUKTBESCHREIBUNG

Colonoscope

Optical System

- Field of View: 120°
- Direction of View: Forward
- Depth of Field: 5 to 100 mm

Insertion Tube

- It Diameter: 13.0 mm

Angulation

- Up/Down: 180°
- Left/Right: 160°



### Lenght

-Working Length: 133 cm.

-Total Length: 160 cm.

### Biopsy Channel

-Channel Size: 3.2 mm

### Compatibility

-Fujinon EPX-201 Video Processor

### Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile

0049 176 205 800 86

Tel / Homezone

0049 40 18 20 79 03

Tel / Büro

0049 40 3 61 61 - 256

Fax.

0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN

**Weight**

10 kg



## GALLERIE

