



## FUJINON EC-590 WI COLONOSCOPE // KOLONOSKOP

**SKU:** S022017

**4.444,00 €**

zzgl. [Versandkosten](#)

Field Of View 140 deg

Depth of Field 3-100 mm

The diameter of the distal end 12.8 mm

The diameter of the flexible part 12.8 mm

Bending Up – 180 deg / down – 180 deg to the  
right – 160 deg / left – 160 deg

Working length 1520 mm

Overall Length 1820 mm

Instrument Channel Diameter 3.8mm instrument  
channel

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Flexible Endoskope & Zubehör](#)



Field Of View 140 deg  
Depth of Field 3-100 mm  
The diameter of the distal end 12.8 mm  
The diameter of the flexible part 12.8 mm  
Bending Up – 180 deg / down – 180 deg to the right – 160 deg / left – 160 deg  
Working length 1520 mm  
Overall Length 1820 mm  
Instrument Channel Diameter 3.8mm instrument channel

Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile	0049 176 205 800 86
Tel / Homezone	0049 40 18 20 79 03
Tel / Büro	0049 40 3 61 61 - 256
Fax.	0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN

<b>Weight</b>	11 kg
---------------	-------



## GALLERIE

