



[GE 6S-D CARDIO-SONDE KARDIOSONDE ECHOSONDE](#)

SKU: AT012024-NV



Categories: [Ultraschall](#) / [Doppler](#) / [Sonden](#) / [Zubehör](#)

PRODUKTBESCHREIBUNG

Array Type: Phased Array

Applikationen: Kinderheilkunde, Kardiologie, Neonatal, Fetales Herz, Abdominal

Frequenzband: 2.4 - 8.0 MHz

Kompatibilität: GE Vivid S60; GE Vivid S70; GE Vivid E9; GE Vivid E90; GE Vivid E95; GE Logiq E9; GE Logiq S7; GE Logiq S8

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

