



[GE 6TV 5.0 MHZ PHASED ARRAY TEE PROBE](#)

SKU: B082014-04

zzgl. [Versandkosten](#)

The GE 6TV ultrasound transducer probe which is compatible with GE Logiq 5, 7, 9 and Vivid 3, 7 ultrasound systems.

It can be used for TEE applications.

2.9 TO 6.7 MHz

14 mm aperture

TEE Array

unter anden für VIVID FIVE

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall / Doppler / Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

The GE 6TV ultrasound transducer probe which is compatible with GE Logiq 5, 7, 9 and Vivid 3, 7 ultrasound



systems.

It can be used for TEE applications.

2.9 TO 6.7 MHz

14 mm aperture

TEE Array

unter anden für VIVID FIVE

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

