



## [GE 6TV 5.0 MHZ PHASED ARRAY TEE PROBE](#)

**SKU:** B082014-04

The GE 6TV ultrasound transducer probe which is compatible with GE Logiq 5, 7, 9 and Vivid 3, 7 ultrasound systems.

It can be used for TEE applications.

2.9 TO 6.7 MHz

14 mm aperture

TEE Array

unter anden für VIVID FIVE

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall](#) / [Doppler](#) / [Sonden](#) / [Zubehör](#)

## **PRODUKTBESCHREIBUNG**

The GE 6TV ultrasound transducer probe which is compatible with GE Logiq 5, 7, 9 and Vivid 3, 7 ultrasound systems.



It can be used for TEE applications.

2.9 TO 6.7 MHz

14 mm aperture

TEE Array

unter anden für VIVID FIVE

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

