



[GE CZB 6.5 MHZ PART NR. 2152424 FOR LOGIQ 200](#)


SKU: I102014-01

CZB Curved Array Transducer 6.5 MHz. For small parts neonatal and biopsy applications.
Part Number 2152424

[Vollfunktionsfähig](#)
[Im guten gebrauchten Zustand](#)
[Artikel-Standort. Hamburg](#)

www.mtc-hamburg.com
info@mtc-hamburg.com

Mobile.	0049 176 205 800 86
Tel / Homezone	0049 40 18 20 79 03
Tel / Büro	0049 40 3 61 61 - 256
Fax.	0049 40 3 61 61 - 257



Categories: [Ultraschall / Doppler / Sonden /Zubehör](#)

PRODUKTBESCHREIBUNG

CZB Curved Array Transducer 6.5 MHz. For small parts neonatal and biopsy applications.
Part Number 2152424

[Vollfunktionsfähig](#)
[Im guten gebrauchten Zustand](#)
[Artikel-Standort. Hamburg](#)



www.mtc-hamburg.com
info@mtc-hamburg.com

Mobile. 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

