



## GE LOGIQ 200 PRO ULTRASOUND

SKU: I102014

1.899,00 €

zzgl. [Versandkosten](#)

GE CAE P9603DB Ultrasound Transducer 5MHz  
GE CZB 6.5 MHz Part Nr. 2152424 for Logiq 200

Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile.	0049 176 205 800 86
Tel / Homezone	0049 40 18 20 79 03
Tel / Büro	0049 40 3 61 61 - 256
Fax.	0049 40 3 61 61 - 257



---

Categories: [Ultraschall](#) / [Doppler](#) / [Sonden](#)  
[/Zubehör](#)

## PRODUKTBESCHREIBUNG

GE CAE P9603DB Ultrasound Transducer 5MHz  
GE CZB 6.5 MHz Part Nr. 2152424 for Logiq 200

Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg



[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

