



[GE RIC 6-12-D ULTRASOUND TRANSDUCER](#)

SKU: Z102018

Brand: GE

Model: RIC 6-12-D

Description: OB/GYN/Urology

Type: Endocavity

Frequency: 5-13MHz

Ultrasound Machine: Voluson E6, Voluson E8

Expert

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall / Doppler / Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

Brand: GE

Model: RIC 6-12-D

Description: OB/GYN/Urology

Type: Endocavity

Frequency: 5-13MHz

Ultrasound Machine: Voluson E6, Voluson E8 Expert



Vollfunktionsfähig
Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

