



## **GE VOLUME 3.5- 11.0 MHZ RSP 5-12 LINEAR VOLUMENSONDE BABYFACE 3D 4D**

**SKU:** 11-032014

Hersteller: GE

Modell : RSP 5-12

Bezeichnung: Volume 3.5- 11.0 MHz

Applikation: Babyface 3D 4D

Für GE Voluson 730 Serie

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall / Doppler / Sonden](#)  
[/Zubehör](#)

## **PRODUKTBESCHREIBUNG**

Hersteller: GE

Modell : RSP 5-12

Bezeichnung: Volume 3.5- 11.0 MHz

Applikation: Babyface 3D 4D

Für GE Voluson 730 Serie

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)



Mobile. 0049 176 205 800 86  
Tel / Homezone 0049 40 18 20 79 03  
Tel / Büro 0049 40 3 61 61 - 256  
Fax. 0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

