



HP 21258B ULTRASOUND TRANSDUCER

SKU: W092014

compatible with HP Sonos 2000 and 2500 ultrasound systems

7.5/5.5MHz

Trap Linear Array

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall](#) / [Doppler](#) / [Sonden](#) / [Zubehör](#)

PRODUKTBESCHREIBUNG

compatible with HP Sonos 2000 and 2500 ultrasound systems

7.5/5.5MHz

Trap Linear Array

Vollfunktionsfähig



Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com

Mobile. 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

