



## [INTERACOUSTICS HEARING AID ANALYZERS](#) [MS25 - AUDIOMETER](#)

**SKU:** B062012

gepflegt  
Vollfunktionsfähig

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)  
0049 40 18 20 79 03  
0049 176 205 800 86  
fax.0049 40 65 86 33 64



---

**Categories:** [Audiometer](#)

## **PRODUKTBESCHREIBUNG**

Features:

- Real Ear Measurements
- Integral Test Box
- 2-CC Reference Coupler

Measurements

- Integral Induction Loop
- Integral 4-colour Printer
- NOAH Compatible
- Dedicated Windows Software available

The MS25 Hearing Aid Analyzer is a fully portable, two-in-one instrument which combines a sophisticated hearing aid test box with a comprehensive real ear measurement system. Using a combination of



dedicated function buttons and software defined "F" keys, easy access to all the test options available on the MS25 is ensured. The MS25 indeed offers a highly

Vollfunktionsfähig  
Im guten gebrauchten Zustand  
Letzte Prüfung 01/2011  
Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)  
0049 40 18 20 79 03  
0049 176 205 800 86  
fax.0049 40 65 86 33 64

specified yet very accessible working platform for both hearing aid testing and hearing aid fitting. An optional keyboard allows relevant patient and hearing aid data to be stored directly in the MS25 and later this can be printed together with relevant test results. Obtained test results can be transferred to an appropriate PC via the built-in RS232C interface. An optional IBM compatible Windows software is available. The MS25 is fully NOAH compatible.

## **ZUSÄTZLICHE INFORMATIONEN**



## GALLERIE

