



## [KAVO K11 KNIEANLASSER TYP 4980 40000 U/MIN](#)

**SKU:** 0042015

**593,81 €**

zzgl. [Versandkosten](#)

KaVo K11 Knieanlasser

Typ 4980

40000 U/min // 200 Watt //

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)


[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Dental und Zahnheilkunde](#)

## **PRODUKTBESCHREIBUNG**

KaVo K11 Knieanlasser

Typ 4980

40000 U/min // 200 Watt //

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

