



MEDISON EC4-9ED ULTRASOUND TRANSDUCER

SKU: AK012014-03

Medison EC4-9/10ED Endo-Cavity Probe
Specification : 4MHz~9MHz / 150° / 10mm

www.mtc-hamburg.com
info@mtc-hamburg.com
0049 40 18 20 79 03
0049 176 205 800 86
fax.0049 40 65 86 33 64



Categories: [Ultraschall / Doppler / Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

Medison EC4-9/10ED Endo-Cavity Probe
Specification : 4MHz~9MHz / 150° / 10mm
Application : OB, Gynecology, Urology

Fits for : Medison SA3200(V1.25.01 or up), SA5000II(V2.12 or up), SA5500, SA6000, SA6000C, SA8800, 128BW, SA8000, PICO, SA9900, X4

Im sehr guten gepflegten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com
0049 40 18 20 79 03
0049 176 205 800 86
fax.0049 40 65 86 33 64



ZUSÄTZLICHE INFORMATIONEN



GALLERIE

