



**MIELE G7736 CD MCU C2 DESINFEKTOR,**  
**THERMODESINFEKTOR, LABOR**  
**SPÜLMASCHINE**

**SKU:** I012015-NV

Laborspülmaschine G 7736 CD, Model C2, MCU-  
Steuerung  
400 V , 03 Phasen, 10,2 Kw,  
Abmessungen 1150 mm x 900 x 700 mm

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile.                   0049 176 205 800 86  
Tel / Homezone   0049 40 18 20 79 03  
Tel / Büro             0049 40 3 61 61 - 256  
Fax.                     0049 40 3 61 61 - 257



---

**Categories:** [Autoklaven und Sterilisatoren](#), [Dental](#)  
[und Zahnheilkunde](#)

## PRODUKTBESCHREIBUNG

Laborspülmaschine G 7736 CD, Model C2, MCU-Steuerung  
400 V , 03 Phasen, 10,2 Kw,  
Abmessungen 1150 mm x 900 x 700 mm

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile.                   0049 176 205 800 86  
Tel / Homezone   0049 40 18 20 79 03  
Tel / Büro             0049 40 3 61 61 - 256



Fax. 0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

