



MIELE G7736 CD MCU C2 DESINFEKTOR,
THERMODESINFEKTOR, LABOR
SPÜLMASCHINE

SKU: I012015-NV

zzgl. [Versandkosten](#)

Laborspülmaschine G 7736 CD, Model C2, MCU-
Steuerung

400 V , 03 Phasen, 10,2 Kw,

Abmessungen 1150 mm x 900 x 700 mm

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Autoklaven und Sterilisatoren](#), [Dental](#)
[und Zahnheilkunde](#)

PRODUKTBESCHREIBUNG

Laborspülmaschine G 7736 CD, Model C2, MCU-Steuerung

400 V , 03 Phasen, 10,2 Kw,

Abmessungen 1150 mm x 900 x 700 mm

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03



Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

