



OLYMPUS BF-1T40 BRONCHOSKOP / FLEXIBLE BRONCHOSCOPE

SKU: CK022016

1.100,00 €

Olympus BF-1T40

Field of view: 120°

Direction of view : 0° (forward viewing)

Depth of view : 3~50mm

Distal End Outer diameter : 5.9mm

Bending Section Range of tip bending : UP 180°, DOWN
130°

Insertion Tube Outer diameter : 6.0mm

Working Length 550mm

Total Length 840mm

Instrument Channel Inner diameter 2.8mm

Minimum visible distance 5mm from distal end

Rückgaberecht . 14 Tage

Versand . Velweit

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257





PRODUKTBESCHREIBUNG

Olympus BF-1T40

Field of view: 120°

Direction of view : 0° (forward viewing)

Depth of view : 3~50mm

Distal End Outer diameter : 5.9mm

Bending Section Range of tip bending : UP 180°, DOWN 130°

Insertion Tube Outer diameter : 6.0mm

Working Length 550mm

Total Length 840mm

Instrument Channel Inner diameter 2.8mm

Minimum visible distance 5mm from distal end

Rückgaberecht . 14 Tage

Versand . Velweit

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

