



OLYMPUS EXERA TJF-160R VIDEO
DUODENOSCOPE VIDEODUODENOSCOPE

SKU: -----

3.200,00 €

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com

| | |
|----------------|-----------------------|
| Mobile | 0049 176 205 800 86 |
| Tel / Homezone | 0049 40 18 20 79 03 |
| Tel / Büro | 0049 40 3 61 61 - 256 |
| Fax. | 0049 40 3 61 61 - 257 |



Categories: [Flexible Endoskope & Zubehör](#)

PRODUKTBESCHREIBUNG

Olympus EXERA TJF-160R Video Duodenoscope.

Removable cap for easy cleaning of the distal end

Superb treatment capability with exceptionally large 4.2mm diameter instrument channel that accommodates a wide range of Endo-Therapy accessories including a variety of stents

Ergonomically designed grip to enhance scope maneuverability and handling

Scope ID function stores individual scope information such as number of procedures performed

Compatible with CV-180, CV-160, CV-140 and CV-100 video processors

Specifications:

Diameter: 11.3mm



Working Length: 124cm

Instrument Channel: 4.2 mm

Field of View: 100 degrees

Angulation Range: Up: 120, Down: 90, Right: 110, Left: 90

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile

0049 176 205 800 86

Tel / Homezone

0049 40 18 20 79 03

Tel / Büro

0049 40 3 61 61 - 256

Fax.

0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

