



## OLYMPUS LF-TP ENDOSKOP INTUBATION SCOPE

**SKU:** AW012016

- \* 1.5mm suction channel
- \* 90° field of view
- \* 3.8mm distal tip and 4.1mm insertion tube
- \* 120° up/down angulation
- \* 600mm working length
- \* Fully immersible for assured reprocessing
- \* Compatible with conventional light sources and video systems for versatility

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Flexible Endoskope & Zubehör](#)

## PRODUKTBESCHREIBUNG

- \* 1.5mm suction channel
- \* 90° field of view
- \* 3.8mm distal tip and 4.1mm insertion tube



- \* 120° up/down angulation
- \* 600mm working length
- \* Fully immersible for assured reprocessing
- \* Compatible with conventional light sources and video systems for versatility

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile.                0049 176 205 800 86  
Tel / Homezone      0049 40 18 20 79 03  
Tel / Büro            0049 40 3 61 61 - 256  
Fax.                    0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

