



## [PENTAX EC-3840F2 VIDEO KOLONOSKOP](#)

SKU: BF012016

**2.499,00 €**

zzgl. [Versandkosten](#)

Field Of View: 120°

IT Diameter : 12.8mm

Angulation Up : 180°

Angulation Down : 180°

Angulation Right : 160°

Angulation Left : 160°

Working Length : 150cm

Channel Size : 3.8mm

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

Categories: [Flexible Endoskope & Zubehör](#)

## PRODUKTBESCHREIBUNG

Field Of View: 120°



IT Diameter : 12.8mm  
Angulation Up : 180°  
Angulation Down : 180°  
Angulation Right : 160°  
Angulation Left : 160°  
Working Length : 150cm  
Channel Size : 3.8mm

Im guten gebrauchten Zustand  
Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)  
[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile.                0049 176 205 800 86  
Tel / Homezone      0049 40 18 20 79 03  
Tel / Büro            0049 40 3 61 61 - 256  
Fax.                    0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

