



PENTAX PVA-1000 VIDEO ENDOSCOPY ADAPTER MIT ADAPTER AL-OF5

SKU: R062015

249,00 €

zzgl. [Versandkosten](#)

PVA-1000 Adapter for use with the Pentax video endoscopy equipment

AL-OF5 Adapter to connect an Olympus OES Endoscope to a Pentax Light Source

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Endoskopie und Zubehör](#)

PRODUKTBESCHREIBUNG

PVA-1000 Adapter for use with the Pentax video endoscopy equipment

AL-OF5 Adapter to connect an Olympus OES Endoscope to a Pentax Light Source

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg



www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

