



## PICKERP/N 800342 3.5MHZ LINEAR PROBE

SKU: P062016

119,00 €

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

Categories: [Ultraschall](#) / [Doppler](#) / [Sonden](#)  
[/Zubehör](#)

## PRODUKTBESCHREIBUNG

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN

**Weight** 2 kg



## GALLERIE

