



## SHIMADZU L123-050D ULTRASOUND TRANSDUCER FOR SHIMADZU SDL-310B

**SKU:** AH102014-02

zzgl. [Versandkosten](#)

Shimadzu L123-050d Linear transducer for shimadzu sld-310b

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall](#) / [Doppler](#) / [Sonden](#)  
[/Zubehör](#)

## PRODUKTBESCHREIBUNG

Shimadzu L123-050d Linear transducer for shimadzu sld-310b

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile

0049 176 205 800 86

Tel / Homezone

0049 40 18 20 79 03

Tel / Büro

0049 40 3 61 61 - 256

Fax.

0049 40 3 61 61 - 257



## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

