



SIEMENS 5MHZ LINEAR SONDE / PROBE

SKU: J032013 & 03-062014 // AJ042014-02

für Sonoline SX/SL/SL1/SL2/CD AC/Si 200/Si250 Si 400 /Si 450

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall / Doppler / Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

für Sonoline SX/SL/SL1/SL2/CD AC/Si 200/Si250 Si 400 /Si 450

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



ZUSÄTZLICHE INFORMATIONEN



GALLERIE

