



SONOACE L5-9 PB-L75/60ED LINEAR ARRAY ULTRASOUND TRANSDUCER PROBE

SKU: F032014-02

zzgl. [Versandkosten](#)

geeignet, unter anderen für .SONOACE // MEDISON
SA-5500; 6000C; 6000-II ; 8000 ...

Frequency 5.0-9 Mhz

Anwendung. Abdomen, Fetal heart, Gynecology and OB

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall / Doppler / Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

geeignet, unter anderen für .SONOACE // MEDISON SA-5500; 6000C; 6000-II ; 8000 ...

Frequency 5.0-9 Mhz

Anwendung. Abdomen, Fetal heart, Gynecology and OB

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03



Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

