



SYSMEX AUTO DILUTOR AD-260

SKU: V042016

999,00 €

zzgl. [Versandkosten](#)

Sysmex AD-260 – Auto Dilutor

Manufactures: Toa, Medical Electronics Co. Ltd.,
Kobe, Japan

Sysmex Auto Dilutor Type AD-260 is used for
diluting blood samples in 2 stages:

1 stage – for WBC determination

2 stage – for RBC determination

Blood is diluted in dilutor with isotonic solution and
next WBC, RBC HCT, HBG and derived parameters
are measured in Cell Counter.

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Laborequipment](#)

PRODUKTBESCHREIBUNG

Sysmex AD-260 – Auto Dilutor

Manufactures: Toa, Medical Electronics Co. Ltd., Kobe, Japan

Sysmex Auto Dilutor Type AD-260 is used for diluting blood samples in 2 stages:

1 stage – for WBC determination

2 stage – for RBC determination

Blood is diluted in dilutor with isotonic solution and next WBC, RBC HCT, HBG and derived parameters are measured in Cell Counter.

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN

Weight

5 kg



GALLERIE

