



TOSHIBA PLF-705S LINEAR ULTRASOUND PROBE/TRANSDUCER

SKU: BT022016 // BZ022016 // CA022016

499,00 €

zzgl. [Versandkosten](#)

Bezeichnung . PLF-705S

Hersteller . TOSHIBA

für Toshiba Tosbee , 140A/250A/270A/340A
systems

Linear Ultraschallkopf

Anwendung : Vaskular

Zustand. Gebraucht //

Artike-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall](#) / [Doppler](#) / [Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

Bezeichnung . PLF-705S



Hersteller . TOSHIBA
für Toshiba Tosbee , 140A/250A/270A/340A systems
Linear Ultraschallkopf
Anwendung : Vaskular

Zustand. Gebraucht //

Artike-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com

Mobile. 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

