



TOSHIBA PLF-705S LINEAR ULTRASOUND PROBE/TRANSDUCER

SKU: AB042013-01 // DH112015-02

Bezeichnung . PLF-705S

Hersteller . TOSHIBA

für Toshiba Tosbee , 140A/250A/270A/340A
systems

Linear Ultraschallkopf

Anwendung : Vaskular

Zustand. Gebrauchte //

Artike-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall / Doppler / Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

Bezeichnung . PLF-705S

Hersteller . TOSHIBA

für Toshiba Tosbee , 140A/250A/270A/340A systems



Linear Ultraschallkopf
Anwendung : Vaskular

Zustand. Gebraucht //

Artike-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com

Mobile. 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

