



**TOSHIBA PLG-805S / PLG 805 S 8 MHZ, 70 MM LINEAR PROBE**

**SKU:** BQ022016

**799,00 €**

zzgl. [Versandkosten](#)

Toshiba Linear 6.0 - 10.0 MHz PLG-805S Schilddrüse  
Mamma Small Part

Toshiba

PLG-805S

Linear 6.0-10.0 MHz

Anwendung : Schilddrüse Mamma Small Part

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



---

**Categories:** [Ultraschall / Doppler / Sonden](#)  
[/Zubehör](#)

## PRODUKTBESCHREIBUNG

Toshiba Linear 6.0 - 10.0 MHz PLG-805S Schilddrüse Mamma Small Part

Toshiba

PLG-805S

Linear 6.0-10.0 MHz



Anwendung : Schilddrüse Mamma Small Part

Artikel-Standort. Hamburg

[www.mtc-hamburg.com](http://www.mtc-hamburg.com)

[info@mtc-hamburg.com](mailto:info@mtc-hamburg.com)

Mobile.           0049 176 205 800 86  
Tel / Homezone   0049 40 18 20 79 03  
Tel / Büro        0049 40 3 61 61 - 256  
Fax.              0049 40 3 61 61 - 257

## ZUSÄTZLICHE INFORMATIONEN



## GALLERIE

