



TOSHIBA PLQ-805A ULTRASOUND TRANSDUCER

SKU: D102015-01 // BZ112017 // A0022018
999,00 €

Toshiba PLQ-805A

Frequency: 6.0-12:0 MHz

Compatible system: Famio 8. ...

Application: Schilddrüse, Mamma, Small Parts ,

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com


info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall / Doppler / Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

Toshiba PLQ-805A

Frequency: 6.0-12:0 MHz

Compatible system: Famio 8. ...

Application: Schilddrüse, Mamma, Small Parts ,



Vollfunktionsfähig
Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

