



**TOSHIBA PVG-366M / PVG 366 M 3,75MHZ,
KONVEX PROBE**

SKU: BW022016 // AS052018

799,00 €

zzgl. [Versandkosten](#)

Toshiba Konvex 3,75 MHz

PVG-366M

Konvex 3.75MHZ

Anwendung : Abdomen

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall](#) / [Doppler](#) / [Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

Toshiba Konvex 3,75 MHz

PVG-366M

Konvex 3.75MHZ

Anwendung : Abdomen

Artikel-Standort. Hamburg



www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

