



[TOSHIBA PVQ-375A ULTRASOUND TRANSDUCER](#)

SKU: BE112017

999,00 €

zzgl. [Versandkosten](#)

Toshiba PVQ-375A 50mm Convex Abdominal
Ultrasound Transducer

Frequency: 3.0-6.0 MHz

Compatible system: Famio 8.

Application: Abdominal

Vollfunktionsfähig

Im guten gebrauchten Zustand

Artikel-Standort. Hamburg

www.mtc-hamburg.com
info@mtc-hamburg.com
Mobile 0049 176 205 800 86
Tel / Homezone 0049 40 18 20 79 03
Tel / Büro 0049 40 3 61 61 - 256
Fax. 0049 40 3 61 61 - 257



Categories: [Ultraschall / Doppler / Sonden](#)
[/Zubehör](#)

PRODUKTBESCHREIBUNG

Toshiba PVQ-375A 50mm Convex Abdominal Ultrasound Transducer

Frequency: 3.0-6.0 MHz

Compatible system: Famio 8.

Application: Abdominal



Vollfunktionsfähig
Im guten gebrauchten Zustand
Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

