



WHALEDENT BIOSONIC UC100 ULTRASONIC CLEANER

SKU: BU112014

399,00 €

zzgl. [Versandkosten](#)

Model UC100

- 115 V - 1.1 A - 60 Hz

- Inside Chamber approx 10.5" X 6" X 2.75"

Vollfunktionsfähig

Artikel-Standort. Hamburg

www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257



Categories: [Dental und Zahnheilkunde](#)

PRODUKTBESCHREIBUNG

Model UC100

- 115 V - 1.1 A - 60 Hz

- Inside Chamber approx 10.5" X 6" X 2.75"

Vollfunktionsfähig

Artikel-Standort. Hamburg



www.mtc-hamburg.com

info@mtc-hamburg.com

Mobile. 0049 176 205 800 86

Tel / Homezone 0049 40 18 20 79 03

Tel / Büro 0049 40 3 61 61 - 256

Fax. 0049 40 3 61 61 - 257

ZUSÄTZLICHE INFORMATIONEN



GALLERIE

